


AFTER 2D GEL

1. Gel-Matching with databases (es: Expasy proteomics Server);
2. Western-Blot with specific antibodies;
3. N-terminal microsequencing;
4. Mass spectrometry: MALDI-TOF, MS-MS


AFTER 2D GEL

1. Gel-Matching with databases (es: Expasy proteomics Server);
2. Western-Blot with specific antibodies;
3. N-terminal microsequencing;
4. Mass spectrometry: MALDI-TOF, MS-MS

Detection via Western Blot



Glioma Silver Stain




**Glioma Western Blot
(anti-p53 antibody)**



1. Gel-Matching with databases (es: Expasy);
2. Western-Blot with specific antibodies;
3. N-terminal microsequencing;
4. Mass spectrometry: MALDI-TOF, MS-MS

Edman degradation

EDMAN DEGRADATION




Phenylisothiocyanate (PTC)
reacts with N-terminal residues



PTH-alanine

A cyclical phenylthiocarbamoyl derivative is formed and hydrolyzed



Peptide shortened by one residue

The thiazolinone amino acid is then selectively extracted and treated with acid to form the more stable phenylthiohydantoin (PTH)- amino acid derivative that can be identified by using chromatography or electrophoresis

Automatic sequencer

Detector




Cartridge 45-48 °C
Camera di conversione 64°C
Colonna cromatografica C18 0.8x250 mm 55°C

RP-HPLC

Sequencing

Protein sequence analysis employing Edman degradation chemistry commonly uses HPLC separation as the means for identification of the PTH-amino acid (phenylthiohydantoin amino acid) produced at each cycle. Each amino acid is red at wavelength 269 nm and is identified by a specific retention time



After 2D gel

1. Gel-Matching with databases (es: Expasy);
2. Western-Blot with specific antibodies;
3. N-terminal microsequencing;
4. Mass spectrometry: MALDI-TOF, MS-MS

Spot identification: peptide mass fingerprint

Spot cutting




Decoloration

“in gel” proteolytic digestion

MS identification

Database searching

YAEKGLQQPVRAYMESEKSEEIDLPK




YAEK	609.31
GLQQPVR	797.45
AYMESEK	857.38
SEEIDLPK	957.54

APOLIPROTEIN A-1 PRECURSOR (APO-A1) (30276.5 Da)

Data	M+H	Delta	start	end	Peptide Sequence	Modifications
submitted	matched	ppm			(Click for Fragment Ions)	
825.4571	825.4583	1.4248	177	183	ROGLVETL(Q)	
1017.5370	1017.5369	0.0782	142	150	QYAPLGEFPR(D)	
1026.5964	1026.5948	1.5881	164	172	QGLPLAQLR(D)	
1040.6121	1040.6104	1.6147	229	237	QKAPYEDLR(Q)	
1218.5796	1218.5755	3.3721	206	217	QKGGGDAKHAH(A)	
1255.6538	1255.6574	2.8906	36	46	QGFATVYVAD(D)	
1260.5928	1260.6013	6.7530	131	139	QVHKEYVPR(D)	
1266.6281	1266.6370	7.0489	120	129	QVPLDIFPR(D)	
1288.6094	1288.6174	6.1758	184	194	ROQLAPSDDLR(Q)	
1305.6310	1305.6439	9.8038	184	194	ROQLAPYEDDLR(Q)	pyroGlu
1388.6975	1388.6963	0.8824	130	139	QKWHKEYVPR(Q)	
1398.6805	1398.6905	7.1643	51	63	RODYVAQFASALGE(Q)	
1482.8101	1482.8208	7.2204	34	46	ROKDAVYVAD(D)	
1576.8138	1576.8213	5.3672	69	82	QKLDNDVILATSLR(V)	
2188.0850	2188.1039	8.6184	83	100	QYREGLGPTVDFWDLR(D)	


Peptide Mass Fingerprinting



Spot recovery

Manual

Keratin contaminations



Automatic

Gel deformation!



Enzymatic digestion

Aim: obtain peptide fragments from the protein of interest

a) **peptide mass fingerprinting**

b) **MS peptide sequencing**


<i>Method of protein cleavage</i>	<i>Site of cleavage</i>	<i>Exception</i>	<i>pH range</i>
Trypsin	C-terminus of R-X, K-X	If X = P	7–9
Endoproteinase Glu-C (V8-DE)	C-terminus of E-X, D-X	If X = P	4–8
Chymotrypsin	C-terminus of F, Y, W, L, I, V, M	If X = P	7.5–8.5
Endoproteinase Lys-C	C-terminus of lysine, K-X	If X = P	8.5–8.8
Arg-C	C-terminus of arginine, R-X	If X = P	7.5–8.5
Elastase	Not very specific. C-terminal side of G, A, S, V, L and I.		8–8.5
Pepsin	C-terminus of F, L and E		2–4
Pronase	Pronase is a mixture of endo- and exo-proteinases. It cleaves almost any peptide bond.		7–8, dependent on proteases present

Complete digestion to obtain the minimum number of missing cleavage sites

Use volatile buffers (for easy elimination after enzymatic digestion)

Why Trypsin?

- Robust, stable enzyme
- Works over a range of pH values & Temp.
- Quite specific and consistent in cleavage
- Cuts frequently to produce “ideal” MW peptides
- Inexpensive, easily available/purified
- Does produce “autolysis” peaks (which can be used in MS calibrations)
 - 1045.56, 1106.03, 1126.03, 1940.94, 2211.10, 2225.12, 2283.18, 2299.18



PEPTIDE MASS FINGERPRINT

Peptide digested with trypsin

YAEKGLQQPVRAYMESEKSEEIDLPK





The 4 derived peptides have a well defined mass

YAEK	609.31
GLQQPVR	797.45
AYMESEK	857.38
SEEIDLPK	957.54

These 4 MW exactly identify the original peptide

After trypsin digestion: peptides are purified by reverse phase chromatography

The high resolving power of Reverse Phase [RP] Chromatography has made it the dominant mode of HPLC for both analytical and preparative separation of peptides and proteins, as well as other biomolecules.





The term “Reverse Phase Chromatography” was used because RP is a form of partition chromatography where **chemically bonded phase** is **hydrophobic or non-polar** (e.g. octadecyl group), and the starting mobile phase (e.g. water) must be more polar than the stationary phase.

This is “reversed” from normal phase chromatography, where the stationary phase is polar or hydrophilic and the starting mobile phase is more non-polar or hydrophobic than the stationary phase, hence the term “Reverse Phase Chromatography”.

After trypsin digestion: peptides are purified with Zip Tip (tips packed with Reversed Phase (C18) resin)



The ZipTip pipette tip is a 10 μL pipette tip with a bed of chromatography media fixed at its end. It is intended for concentrating and purifying peptide, protein or oligonucleotide samples. C_{18} and C_4 reversed-phase media for desalting and concentrating peptides and proteins.



How to detect peptides?

After digestion sample is ionized by:


- Protonation $M + H^+ \rightleftharpoons MH^+$
- Cationization $M + Cat^+ \rightleftharpoons M Cat^+$
- Deprotonation $MH \rightleftharpoons M^- + H^+$
- Electron release $M \rightleftharpoons M^+ + e^-$
- Electron capture $M + e^- \rightleftharpoons M^-$

Analytes are ionized in gaseous phase

~~“A MASS SPECTROMETER
MEASURES THE MW ...”~~


“...A MS ANALYSIS GIVES
THE MASS-TO-CHARGE RATIO (m/z)
OF IONS...IN GAS PHASE”.

Mass spectrometer scheme




Mass spectrum: intensity vs m/z ratio plot representing a chemical analysis

Mass spectrometer scheme




Typical Mass Spectrum

- Characterized by sharp, narrow peaks
- X-axis position indicates the m/z ratio of a given ion (for singly charged ions this corresponds to the mass of the ion)
- Height of peak indicates the relative abundance of a given ion (not reliable for quantitation)
- Peak intensity indicates the ion's ability to desorb or “fly”






Maldi ToF mass spectrum of the tryptic digest of a high molecular weight (85,000) glutenin (DY10)


Mass Spectrometer Schematic



Different Ionization Methods

- Electron Impact (**EI** - Hard method)
 - small molecules, 1-1000 Daltons, structure
 - Fast Atom Bombardment (**FAB** – Semi-hard)
 - peptides, sugars, up to 6000 Daltons
 - Electrospray Ionization (**ESI** - Soft)  1980's
 - peptides, proteins, up to 200,000 Daltons
 - Matrix Assisted Laser Desorption (**MALDI**-Soft)  1985
 - peptides, proteins, DNA, up to 500 kD
- 
- ✓ **Soft ionization techniques keep the molecule of interest fully intact**
 - ✓ Made it possible to analyze large molecules via inexpensive mass analyzers such as quadrupole, ion trap and Time-of-Flight (TOF)

FAB - Fast Atom Bombardment




The sample is mixed with a matrix (non-volatile chemical protection environment) and is bombarded under vacuum with a high energy (4000 to 10,000 electron volts) beam of atoms. The atoms are typically from an inert gas such as argon or xenon. Common matrices include glycerol, thioglycerol.

FAB is a relatively soft ionization technique and produces primarily intact protonated molecules denoted as $[M+H]^+$ and deprotonated molecules such as $[M-H]^-$. **Suitable for peptides, small proteins (up to 5000 Da)**


Soft Ionization Methods

337 nm UV laser



MALDI

Matrix Assisted Laser Desorption




ESI

Electrospray Ionization


ESI

The ESI source operates at atmospheric pressure. A sample solution is sprayed from a small tube into a strong electric field in the presence of a flow of warm nitrogen to assist desolvation. The droplets formed evaporate in a region maintained at a vacuum of several torr causing the **charge to increase on the droplets**. The multiply charged ions then enter the analyzer.



ions carry multiple charges, which reduces their mass-to-charge ratio compared to a singly charged species. This allows mass spectra to be obtained for large molecules

ESI suitable for combination with HPLC and CE (Capillary Electrophoresis)




When the charge density at the droplet surface reaches a critical value (the Rayleigh limit), a so called **Coulombic explosion** occurs and several even smaller droplets are formed, each carrying some fraction of the original droplets surface charge.

Solvent evaporation
Droplet contraction
Coulombic explosions

repeated until the molecular adducts are released from the final droplet.

If a positive voltage is applied to the capillary, then the droplets will carry positive charges and finally positive ions are formed, such as $[M+H]^+$ and $[M+Na]^+$ adducts. In the negative-ion mode, the base peak is typically the $[M-H]^-$ ion.



W.M. Keck FT-ICR Mass Spectrometry Laboratory

Electrospray Ionization


- Can be coupled to HPLC
- Can be modified to “nanospray” system with flow rates $< 1 \mu\text{L}/\text{min}$
 - normal ESI flow rates $1\text{-}500 \mu\text{L}/\text{min}$
 - normal analytical HPLC $0.1\text{-}1 \text{ ml}/\text{min}$
 - nanoESI $1\text{-}500 \text{ nl}/\text{min}$
- Very sensitive technique, requires less than a picomole of material
- Strongly affected by salts & detergents
- Positive ion mode measures $(M + H)^+$ (add formic acid to solvent)
- Negative ion mode measures $(M - H)^-$ (add ammonia to solvent)

Positive or Negative Ion Mode?

- If the sample has functional groups that readily **accept H⁺** (such as amide and amino groups found in **peptides and proteins**) then **positive ion** detection is used –
For positive-ion mode, 0.1% formic acid or acetic acid is usually added into the analyte solution to enhance protonation and increase sensitivity-
- If a sample has functional groups that readily **lose a proton** (such as carboxylic acids and hydroxyls as found in **nucleic acids and sugars**) then **negative ion** detection is used
For negative-ion mode, 0.3% NH₄OH is usually added into the analyte solution to help deprotonation and increase sensitivity -


Multiply Charged Ions

- Samples of MW up to 1200 Da usually produce singly charged ions (parent mass + H)
- Larger samples (typically peptides) yield ions with multiple charges (from 2 to 20 +)
- Multiply charged species form a Gaussian distribution with those having the most charges showing up at lower m/z values



ESI spectrum of Lysozyme
MW = 14,305.14

➡
deconvolution




ESI Transformation

- Software can be used to convert these multiplet spectra into single (zero charge) profiles which gives MW directly
- This makes MS interpretation much easier and it greatly increases signal to noise

Two methods are available

- Transformation (requires prior peak ID)
- Maximum Entropy (no peak ID required)


Maximum Entropy




ESI and Protein Structure

- ESI spectra are actually quite sensitive to the **conformation** of the protein
- Folded, ligated or complexed proteins tend to display **non-gaussian peak distributions**, with few observable peaks weighted toward higher m/z values
- Denatured or open form proteins/peptides which ionize easier tend to display many peaks with a classic gaussian distribution


Native Azurin



Denatured Azurin




MALDI Ionization




- Absorption of UV radiation by chromophoric matrix and ionization of matrix
- Dissociation of matrix, phase change to super-compressed gas, charge transfer to analyte molecule
- Expansion of matrix at supersonic velocity, analyte trapped in expanding matrix plume (explosion/"popping")

MALDI


- Unlike ESI, MALDI generates spectra that have just a **singly charged ion**
- Positive mode generates ions of $M + H$
- Negative mode generates ions of $M - H$
- Generally **more robust than ESI** (tolerates salts and nonvolatile components)
- Easier to use and maintain, capable of higher throughput
- Requires 10 μL of 1 $\text{pmol}/\mu\text{L}$ sample





MALDI Sample Limits

- Phosphate buffer < 50 mM
- Ammonium bicarbonate < 30 mM
- Tris buffer < 100 mM
- Guanidine (chloride, sulfate) < 1 M
- Triton < 0.1%
- SDS < 0.01%
- Alkali metal salts < 1 M
- Glycerol < 1%



....MALDI or Electrospray ?

ESI is **better for the analysis of complex mixture** as it is directly interfaced to a separation technique (i.e. HPLC or CE)

MALDI is more “flexible” (MW from 200 to 400,000 Da)

MALDI is easier to use and maintain and more robust (allows analysis of bacteria, cells and tissues without purification)

MALDI

ESI

<u>Sensitivity</u> =	femtomole 10^{-15} M (...attomole 10^{-18} M)	
<u>Simplicity</u> =	very easy	training required
<u>\$\$\$</u> =	70 to 650 k\$	120 to 650 k\$
<u>Speed</u> (high throughput)=	$\sim 10^4/\text{day}$	dynamic system

MASS ACCURACY AND RESOLUTION

Mass accuracy in mass spectrometry is related to the **calibration of the mass analyzer** to properly assign the true mass-to-charge ratio to a detected ion and to the **resolution of the detector response**, which is in the form of an intensity spike or peak within a mass spectrum. In calibrating a mass spectrometer, a series of **standard compounds**, usually in the form of a multicomponent standard, are measured and related to the fundamental properties of the mass analyzer.

The mass accuracy in mass spectrometry is usually calculated as a parts-per-million (ppm) error where the theoretical mass (calculated as the monoisotopic mass*) is subtracted from the observed mass, divided by the observed mass and multiplied by a 10^6 factor. Mass accuracy (ppm error) can thus be represented as

$$\text{Mass accuracy (ppm)} = \frac{m_{\text{observed}} - m_{\text{theoretical}}}{m_{\text{observed}}} \times 10^6$$

*Monoisotopic Mass: the sum of the exact or accurate masses of the lightest stable isotope of the atoms in a molecule

Mass error

- Mass error = (exact mass) – (accurate mass)
- Mass error in parts per million (ppm) =


$$\frac{\text{(mass error)}}{\text{(exact mass)}} \times 10^6$$

MASS RESOLUTION

In mass spectrometry, mass resolution is generally calculated as $m/\Delta m$ where m is the m/z value obtained from the spectrum and Δm is the full peak width at half maximum (FWHM).

In the figure for the m/z 703.469 ion, which has a full peak width at half maximum of 0.1797 m/z , the resolution of the m/z 703.469 would be calculated as $m/\Delta m = 703.469/0.1797 = 3914$.

The narrower and sharper a peak is will result in a higher-resolution value and a better estimate of the apex of the Gaussian-shaped peak, which thus results in a better estimate of the true value of the mass-to-charge ratio



Full width at half maximum (FWHM) calculation of the resolution of m/z 703.469 as $m/\Delta m = 3914$.

Mass resolution


Mass resolution = (ion mass)/(mass peak width)

$$R = (m/z) / W_{1/2}$$

$W_{1/2}$ = Peak width
at half height
 $W_{1/2} = 0.5 \text{ Da.}$

Example:
 $R = 200 \text{ Da.} / 0.5 \text{ Da.}$
 $R = 400$

$m/z = 200 \text{ Da.}$






FIGURE 28.15 Effect of mass resolving power on the mass spectrum.

A mass resolution of approximately 200 for the overall curve where none of the isotopic peaks have been resolved. A resolution of 200 will give an average mass value for the species being measured with a high degree of error for the exact mass (>500 ppm).

The three isotope peaks are resolved at a resolution of 2000, and an intermediate accurate exact mass is obtainable.


At a resolution of 25,000, the three isotopic peaks are completely baseline resolved, and a mass accuracy typically less than 5 ppm can be obtained

MALDI resolution issues

There are two known problems with MALDI ionization in the gaseous plume of analytes above the target that cause a decrease in the resolution of the detected analytes

1. initial spatial distribution where not all of the desorbed analytes are at the exact same distance from the detector at the start of their flight (placement in crystal structure may attribute to this).
2. Not all of the analytes may have exactly the same velocity at the start of their flight toward the detector (*see slides on TOF, REFLECTRON, improvements from $m/\Delta m = 500$ to $m/\Delta m = 10,000$ or higher*).

Mass Spectrometer Schematic



Mass analyzers

The mass analyzer is the heart of the mass spectrometric instrumentation used in the separation of molecular ions (M^+) and analyte ions (e.g., $[M+H]^+$) in the gas phase.


The most fundamental aspect of the mass analyzer is the ability to separate ions according to their mass-to-charge (m/z) ratio

Mass Analyzers

- Magnetic Sector Analyzer (MSA)
 - High resolution, exact mass, original MA
- Quadrupole Analyzer (Q)
 - Low resolution, fast, cheap
- Time-of-Flight Analyzer (TOF)
 - No upper m/z limit, high throughput
- ORBITRAP
 - high sensitivity, resolution and mass accuracy
- Ion Trap Mass Analyzer (QSTAR)
 - Good resolution, all-in-one mass analyzer
- Ion Cyclotron Resonance (FT-ICR)
 - Highest resolution, exact mass, costly

Magnetic sector mass analyzers

When charged particles enter a magnetic field, they will possess a circular orbit that is perpendicular to the poles of the magnet




The ions enter a flight tube (first field-free region) from a source through a source exit slit and travel into the magnetic field. The accelerating voltage in the source will determine the kinetic energy (KE) that is imparted to the ions:

$$KE = zeV = \frac{1}{2}mv^2$$

V is the accelerating voltage in the source, e is the fundamental charge of an electron ($1.60 \times 10^{-19} \text{ C}$), m is the mass of the ion, v is the velocity of the ion, and z is the number of charges. The magnetic field will deflect the charged particles according to the radius of curvature of the flight path (r) that is directly proportional to m/z of the ion

Quadrupole Mass Analyzer

- A quadrupole mass filter consists of four parallel metal rods with different charges
- Two opposite rods have an applied potential of $(U+V\cos(wt))$ and the other two rods have a potential of $-(U+V\cos(wt))$
- The applied voltages affect the trajectory of ions traveling down the flight path
- For given dc and ac voltages, **only ions of a certain mass-to-charge ratio pass through the quadrupole filter** and all other ions are thrown out of their original path







FIGURE 30.11 Quadrupole orientation and the configuration for the connections of the DC voltage (U) and radio frequency (rf) voltage (V). Ions are accelerated into the quadrupole by a small voltage of 5 eV, and under the influence of the combination of electric fields, the ions follow a complicated trajectory path.

Quadrupole Mass Analyzer



Quadrupole Mass Analyzer

Ions are accelerated into the quadrupole by a small voltage of 5 eV, and under the influence of the combination of electric fields, the ions follow a complicated trajectory path. If the oscillation of the ions in the quadrupole has finite amplitude, it will be stable and pass through. If the oscillations are infinite, they will be unstable and the ion will collide with the rods.

The DC and rf voltages have been selected to give an m/z value of 100 a stable trajectory through the quadrupoles. An m/z value of 10, which is a less massive ion, will have a very unstable trajectory and will collide with the quadrupole rods at an early stage. An m/z 1000 species will tend to travel further through the quadrupole field but still possesses an unstable trajectory and will also suffer collision with the rods and thus also be effectively filtered out.




FIGURE 30.12 Stable and unstable trajectories of ions through the quadrupole. The m/z 100 species has been selected for stable path and transmission through the quadrupole for detection.

ORBITRAP


Ions are trapped in an electrostatic field between an inner and outer electrode (Makarov, 2000).

As the ions rotate around the inner electrode, they precess along its axis with a frequency characteristic of their mass-tocharge (m/z) ratio.

Acquisition of transients and the Fourier transformation of that signal yields frequencies and their intensities.

A simple relationship converts frequencies into m/z values.


It features a pair of axially symmetric electrodes: a central spindle-like electrode and an outer «barrel-like electrode». In this electric field ions rotate around the central electrode while oscillating down the length of the electrode. The frequency of such oscillations is proportional to $(m/z)^{-1/2}$



TOF

The most common mass spectrometer that is coupled to the MALDI ionization technique is the time-of-flight (TOF) mass spectrometer (TOF/MS). The TOF mass spectrometer separates compounds according to their mass-to-charge (m/z) ratios through a direct relationship between a compound's drift time through a predetermined drift path length and the analyte ion's mass-to-charge (m/z) ratio.

Principle of Time-of-Flight (TOF)



Initially, all the ions have similar KE imparted to them from the draw-out pulse (representing time zero), which accelerates them into the flight tube. Because the compounds have different masses, their velocities will be different according to the relationship between KE and mass represented by $KE = zeV = 1/2mv^2$. From this expression, the mass-to-charge ratio is related to the ion's flight time by the following expression: $m/z = 2eVt^2/L^2$:

$$KE = zeV = \frac{1}{2}mv^2$$

$$v = \left(\frac{2zeV}{m} \right)^{\frac{1}{2}}$$

$$t = \frac{L}{v}, \quad L = \text{length of drift tube}$$

$$t = L \left(\frac{m}{2zeV} \right)^{\frac{1}{2}}, \quad V \text{ and } L \text{ are fixed}$$

Solving for m/z ,

$$t^2 = L^2 \left(\frac{m}{2zeV} \right)$$

$$\frac{m}{z} = \frac{2eVt^2}{L^2}$$

Mass Spec Equation (TOF)

$$\frac{m}{z} = \frac{2Vt^2}{L^2}$$

m = mass of ion

z = charge of ion


V = voltage

L = drift tube length


t = time of travel

Normal and Reflectron analyzers

Analizzatore a tempo di volo lineare



Analizzatore a tempo di volo a riflettometro





The mass spec resolution can be improved by increasing the length of the flight tube or (better) by refocusing the ion beam with a “reflectron”

Ions with same m/z ratio but different kinetics energies are refocused by the reflectron therefore reducing differences in flight time due to different initial velocities

Reflectron

M@LDI™ LR by Micromass, UK






FIGURE 30.7 Components of a time-of-flight mass spectrometer illustrating the major sections including the source, drift tube, reflectron electrostatic mirror, and detector.







FIGURE 30.9 Electrostatic mirror focusing of two ions that have the same m/z value but slightly different kinetic energies. Ion path 1 possesses slightly higher kinetic energy in relation to ion path 2. Ion path 1 travels slightly farther to match that of ion path 2. The two ions are focused and arrive at the detector at the same time.

Mass Spectrometer Schematic



MS Detectors

- Early detectors used photographic film
- Today's detectors (**ion channel and electron multipliers**) produce electronic signals via secondary electronic emission when struck by ion
- Timing mechanisms integrate these signals w scanning voltages to allow the instrument to report which m/z has struck the detector
- Need constant and regular calibration



Electron Multiplier (Dynode)

Hybrid instruments /Tandem MS

Combines two or more mass analyzers of the same or different types

First mass analyzer isolates the ion of interest (parent ion)

The ions are then fragmented between the first and second mass analyzer via collisions

The last mass analyzer obtains the mass spectrum of the fragments ions (daughter ions spectrum)


MS-MS spectra reveal fragmentation patterns to provide structural information about a molecule

- Different MS-MS configurations
 - Quadrupole-quadrupole (low energy)
 - Magnetic sector-quadrupole (high)
 - Quadrupole-time-of-flight (low energy)
 - Time-of-flight-time-of-flight (low energy)


Tandem Mass Spectrometry

- Purpose is to fragment ions from parent ion to provide structural information about a molecule
- Also allows separation and identification of compounds in complex mixtures
- Uses two or more mass analyzers/filters separated by a collision cell filled with Argon or Xenon
- Collision cell is where selected ions are sent for further fragmentation

Triple quadrupole mass spectrometer




Q-TOF Mass Analyzer



QTOF has 3 main advantages :


- Better mass accuracy (ID proteins/peptides)
- Better resolution (very useful for multiply charged species from ESI)
- Better sensitivity for product ion “scan”, due to duty cycle effects

CID: Collision-Induced Dissociation
CAD: Collision Activated Dissociation




Ion enters into collision cell

CID: Collision-Induced Dissociation
CAD: Collision Activated Dissociation




CID: Collision-Induced Dissociation

CAD: Collision Activated Dissociation



CID: Collision-Induced Dissociation


CAD: Collision Activated Dissociation



Fragmentation of the ion (to decrease E)


CID: Collision-Induced Dissociation

CAD: Collision Activated Dissociation



Ions produced by fragmentation can be analyzed by second analyzer


ESI-QTOF MS



Proteomics Applications

- Protein sample identification/confirmation
- Protein sample purity determination
- Detection of post-translational modifications
- Detection of amino acid substitutions
- De novo peptide sequencing
- Mass fingerprint identification of proteins
- Monitoring protein folding
- Monitoring protein-ligand complexes/struct.

Bottom-up approach



In the **bottom-up approach**, proteins in complex mixtures can be **separated before enzymatic (or chemical) digestion** followed by direct peptide mass fingerprinting-based acquisition or further peptide separation on-line coupled to tandem mass spectrometry.

Alternatively, the protein mixture can be directly digested into a collection of **peptides** ('shotgun' approach), which are then separated by multidimensional chromatography on-line coupled to tandem mass spectrometric analysis.

Digestion

Peptide Mixture

2D Chromatography


RP

SCX

> 1,000 Proteins Identified

SEQUEST® DTASelect & Contrast

MS/MS Spectrum



Software for protein identification


SEQUEST was one of the first tandem mass spectrometry database search program used for protein identification that correlates uninterpreted tandem mass spectra of peptides with amino acid sequences from protein and nucleotide databases.

MASCOT <http://www.matrixscience.com/>

Proteome Discoverer Software

Multiple database search capability provides the option of applying multiple search algorithms (SEQUEST, Z-Core, Mascot, etc.) and combining their outputs to maximize and cross-validate results

Top-down approach




In the **top-down approach**, proteins in complex mixtures are fractionated and separated into pure single proteins or less complex protein mixtures, followed by off-line static infusion of sample into the mass spectrometer for **intact protein mass measurement** and **intact protein fragmentation**. An on-line LC–MS strategy can also be used for large-scale protein interrogation.

Software for protein identification

- Mascot is a powerful search engine which uses mass spectrometry data to identify proteins from primary sequence databases.
- X! Tandem open source is software that can match tandem mass spectra with peptide sequences

Peptide Mass Fingerprinting by MALDI-TOF...

...the most common procedure to identify proteins by MS



MS analysis of proteolytic fragments is a common way to identify a protein.

The following masses are entered for protein identification: 737.99 - 874.44 - 936.48 - 1030.65 - 1047.06 - 1153.67 - 1269.76 - 1428.85 1536.76 - 1676.97 - 1808.06 - 1862.99 - 2163.26 - 2274.28


Low MW peaks (e.g. < 500 Da) are not generally used because of the high interference of the matrix.

Peptide Mass Fingerprinting

- Used to identify protein spots on gels or protein peaks from an HPLC run
- Depends on the fact that if a peptide is cut up or fragmented in a known way, the resulting fragments (and resulting masses) are unique enough to identify the protein
- Requires a database of known sequences
- Uses software to compare observed masses with masses calculated from database

Peptide Mass Fingerprint (PMF)

Database search



Match between theoretical mass and experimental ones

MASCOT

<http://www.matrixscience.com>

ALDENTE

<http://www.expasy.org/tools/aldente/>

PROFOUND


<http://www.unb.br/cbsp/paginiciais/profound.htm>

MS PROSPECTOR

<http://prospector.ucsf.edu/prospector/4.0.7/html/msfit.htm>

SWISS-PROT

<http://www.expasy.org/sprot/>



Mass list

1428.482505
1434.487611
1675.51502
1803.562355
2123.519746
2336.502887
2464.526762

ALDENTE: sample

All **Sample** Protein Peptide Thresholds Scoring Output Settings About No help

SAMPLE

Sample name unknown

Mw

pl

Peak list

874.479196
993.45726
1007.475247
1045.614758
1428.718073
1631.332142
1675.872109
1676.874966
1753.920078

Upload a file

File type

Mass list

Clear Remove known masses View spectrum

Submit

Example: This can be obtained from a 2-D gel.

Each peak is a pair (mass, intensity). Ideally each peak should represent a peptide in the sample. This is not the case because some peaks are not present, some are noise, besides chemical or physical modifications can occur on peptides changing their masses. Often proteins are isolated using 2D gels, this gives an idea of the pl and Mw of the unknown proteins. If known, you can specify Mw and pl for your sample.

ALDENTE: protein

All	Sample	Protein	Peptide	Thresholds	Scoring	Output	Settings	About	No help
-----	--------	----------------	---------	------------	---------	--------	----------	-------	---------

Database(s)	UniProtKB/Swiss-Prot UniProtKB/TrEMBL Contaminants
Predefined Taxon	All
NCBI TaxID(s)	
Mw min	10000
Mw max	20000
pI min	0
pI max	14
Varsplices	<input checked="" type="checkbox"/>
Fragments	<input checked="" type="checkbox"/>
Keyword(s)	
Keywords mode	Or
Limit to AC	

Define specific NCBI TaxID(s) to search for. Use [Newt](#) database to

Submit

SwissProt is a protein sequence database which provides a high level of annotation (such as the description of the function of a protein, its domains structure, post-translational modifications, variants, etc.), a minimal level of redundancy and high level of integration with other databases

TrEMBL (Traslation from EMBL) is a supplement of Swiss-Prot that contains all the translations of EMBL nucleotide sequence entries not yet integrated in Swiss-Prot.

ALDENTE: peptide

All	Sample	Protein	Peptide	Thresholds	Scoring	Output	Settings	About	No help
-----	--------	---------	----------------	------------	---------	--------	----------	-------	---------

Enzyme

Missed cleavage

Resolution

Ion mode

Digestion

Trypsin

1

Monoisotopic

[M+H]

PTMs

Use PTM

☒

Modifications

Carboxyamidomethyl

Add to list

Remove	Label	Locus	Add	Supp	Mode	Max
<input checked="" type="checkbox"/>	CAM	C	C2H4ON	H	Fixed	1
<input type="checkbox"/>	MSO	M	O		Variable	2

Submit

THEORETICAL PEPTIDES

Enzyme:
Specify the enzyme (see [list and cleavage rules](#)) you used to generate your peptides.

Missed cleavage:
Select the number of missed cleavages allowed.

Resolution:
Specify the isotopic resolution of the experimental masses.
The theoretical masses of the peptides will be calculated accordingly.

Ion mode:
Specify the charge state of the peptides:
- Protonated molecular ions,
[M+H]⁺


* CAM carboxyamidomethyl cysteine (produced by the reaction of **iodoacetamide** with cysteine groups)

ALDENTE: thresholds

All Sample Protein Peptide **Thresholds** Scoring Output Settings About No help

$$\text{ppm} = 10^6 \times \Delta m / m$$

Δm = mass tolerance in ppm



Spectrometer shift max	±	<input type="text" value="1"/>	(Dalton)
Spectrometer slope max	±	<input type="text" value="200"/>	(ppm)
Spectrometer internal error max	±	<input type="text" value="25"/>	(ppm)
Minimum number of Hits		<input type="text" value="4"/>	

To permit detection of mass peaks even in the presence of a slightly erroneous calibration

THRESHOLDS

The calibration of the spectrometer is assimilated to an affine function (a line). You have to define thresholds for this calibration : shift and slope max.

Shift max:

Defines the maximum difference that you allow between a theoretical peptide and an experimental peak in Daltons (absolute value).

Slope max:

Defines the maximum difference that you allow between a theoretical peptide and an experimental peak in ppm (relative value "parts per million").

You can set Shift or Slope to zero.

Submit

ALDENTE: scoring

All	Sample	Protein	Peptide	Thresholds	Scoring	Output	Settings	About	No help
-----	--------	---------	---------	------------	----------------	--------	----------	-------	---------

Peptide Scoring

Intensity	<input type="text" value="1"/>	<input checked="" type="checkbox"/>
Missed cleavage	<input type="text" value="0.5"/>	<input checked="" type="checkbox"/>
C-term R	<input type="text" value="1.5"/>	<input checked="" type="checkbox"/>
C-term K	<input type="text" value="0.5"/>	<input checked="" type="checkbox"/>
C-term limit	<input type="text" value="5"/>	<input checked="" type="checkbox"/>
Modif "CAM"	<input type="text" value="0.3"/>	<input checked="" type="checkbox"/>
Modif "MSO"	<input type="text" value="0.9"/>	<input checked="" type="checkbox"/>

Protein Scoring

Coverage	<input type="text" value="1"/>	<input checked="" type="checkbox"/>
Mw	<input type="text" value="1"/>	<input checked="" type="checkbox"/>
pI	<input type="text" value="1"/>	<input checked="" type="checkbox"/>

Statistics on random Sequences ☒

SCORING

You can check or uncheck the different parameters to be used in the score.

There are two levels in the score :

a) scoring all theoretical peptides and b) scoring the global protein.

Final score = Scoring peptides * Scoring protein

Scoring peptides = Sum of each peptide score. (Default peptide score = 1)

Scoring protein = Scoring Coverage * Scoring pI * Scoring Mw

Intensity:
The peptide score will be multiplied

ALDENTE: output

All	Sample	Protein	Peptide	Thresholds	Scoring	Output	Settings	About	No help
-----	--------	---------	---------	------------	---------	---------------	----------	-------	---------

Maximum number of proteins to display

Merge protein with sequence similarity greater than

Output format

Sort peptides (HTML format)

Your E-mail address

OUTPUT

Number of results:
Define the maximum number of proteins to be displayed.

Format:
Several output formats are available.

Sort peptides (HTML format):
For the HTML format, you can sort peptides either by position on the sequence or by mass values.

E-mail:
Results can be sent by e-mail. Enter your E-mail address in the text area. If no e-mail is provided, the results are displayed on-line in your browser window.
Note! Your E-mail server should allow you to receive HTML.

ALDENTE: results

The results are composed of three parts:

The Header gives general information on the request

The Summary displays a summary of the best proteins candidates

The details are displayed for each protein candidate

Sample	unknown / Peaks 9 / Mass [874.479196 - 1753.920078] / Intensity [1 - 1]
Date	29/11/2006 09:35:09 UTC
Release	UniProtKB/Swiss-Prot Release 51.1 of 14-Nov-2006: 251512 entries.
Proteins	- In range 46923 / After digestion 40 - First Analysis on 40 sequences : After Alignment 1 - Second Analysis on best 1 of first analysis : After Alignment 1 / Displayed 1
Peptides	Generated 1901717 / Matching a peak 13873 / Average per protein 40
Random	Generated 46923 / Best score 0.40 / Mean score 0.00 / Standard deviation 0.00

In range: Number of sequences in the selected database(s) and taxon(s), in the defined mass range, pI range, keywords, etc... Several sequences can be merged in a single Swiss-Prot or TrEMBL entry, whereas sequence variants for one protein entry are counted as separated sequences.

After digestion: Number of proteins with at least the number of hits after digestion and before the alignment process. This means: Number of proteins with enough theoretical peptides matching an experimental peak in the error space defined by the user to limit the comparison (e.g. 0.2 Da and 200 ppm) and before alignment (e.g. 25 ppm for internal error).

First analysis: The first analysis is run with a low precision for a quick scan of the set of proteins selected according to the previous filters (database, taxon, mass and pI ranges, number of hits, error tolerance). This first run will discarded very unlikely proteins.

After alignment: Number of proteins with at least the number of hits aligned (low resolution).

Second analysis: The second analysis is run on the best X proteins of the first analysis (X is the number of proteins to displayed, at least 30). If the number of proteins after the first analysis is less than 30, only these proteins goes into the second analysis.

After alignment: Number of proteins with at least the number of hits aligned (high resolution).

Displayed: Number of proteins displayed defined by the user or less if there are not enough proteins to display.

Note! Because of the higher resolution, some proteins can be discarded after the second analysis. Then the protein number displayed may be less than the max number of proteins defined by the user.

Summary

Identification Summary

Rank	Score	Hits	AC	ID	DE	Mw	pI	Cov	TaxId
					30 first characters	kDa		%	
<u>1</u>	65.75	21	P17844	DDX5_HUMAN (C_1)	Probable ATP-dependent RNA hel...	69	9.1	33	9606
<u>2</u>	32.37	13	Q9Y2X3	NOL5_HUMAN (C_1)	Nucleolar protein 5.	60	9.0	28	9606
<u>3</u>	7.27	6	P00761	TRYP_PIG (C_1)	Trypsin.	23	8.3	31	9823
<u>4</u>	3.74	9	P02769	ALBU_BOVIN (C_1)	Serum albumin.	66	5.6	18	9913
<u>5</u>	2.53	10	Q8K4G5	ABLM1_MOUSE (C_1)	Actin-binding LIM protein 1.	97	8.9	13	10090
<u>6</u>	1.92	6	P55194	3BP1_MOUSE (C_1)	SH3 domain-binding protein 1.	65	5.5	14	10090
<u>7</u>	1.87	10	Q8K449	ABCA9_MOUSE (C_1)	ATP-binding cassette sub-famil...	183	6.5	7	10090
<u>8</u>	1.83	13	P41233	ABCA1_MOUSE (C_1)	ATP-binding cassette sub-famil...	254	6.4	7	10090
<u>9</u>	1.63	6	Q9Z176	2A5R_MOUSE (C_1)	Protein phosphatase 2A, 59 kDa...	56	4.9	13	10090
<u>10</u>	1.61	7	Q6P542	ABCF1_MOUSE (C_1)	ATP-binding cassette sub-famil...	95	6.2	11	10090
<u>11</u>	1.59	6	P12023	A4_MOUSE (C_1)	Amyloid beta A4 protein.	85	4.7	11	10090
<u>12</u>	1.55	9	Q8K442	ABC8A_MOUSE (C_1)	ATP-binding cassette sub-famil...	184	6.9	6	10090
<u>13</u>	1.53	8	Q80XL6	ACD11_MOUSE (C_1)	Acyl-CoA dehydrogenase family ...	87	8.7	11	10090
<u>14</u>	1.39	7	Q99P81	ABCG3_MOUSE (C_1)	ATP-binding cassette sub-famil...	74	6.7	13	10090
<u>15</u>	1.27	6	P48193	41_MOUSE (C_1)	Protein 4.1.	96	5.4	8	10090

Score

Aldente score of the protein. The background color indicates in red or green if the score is lower or greater than the best random score; this threshold should help to eliminate wrong identifications, taking into account the possibility to produce a better score randomly for these proteins.

- Details of a candidate

1) [P00698](#) LYSC_CHICK (C_1) [Up](#)

Lysozyme C.

UniProtKB/Swiss-Prot: [9031](#)

Score	Z-score	Mw	pI	Hits	Coverage	Shift (Da)	Slope (ppm)
18.42	9872.50	14313	9.33	6	51%	0.04	17

methionine oxidation (**MSO**)
carbamidomethyl cysteine (**CAM**)

Exp	Theo	Intensity			Delta	Dev	MC	CAM	MSO	PTM	Position		Sequence
Da	Da	UI	%	rank	Da	ppm	ppm				start	end	
874.479196	874.416572	1	100	1	0.06	72	9	-	-	-	33	- 39	HGLDNYR
1428.718073	1428.6502	1	100	1	0.07	48	3	-	-	-	52	- 63	FESNFNTQATNR
1753.920078	1753.835092	1	100	1	0.08	48	9	-	-	-	64	- 79	NTDGSTDYGILQINSR
* 993.45726	993.39954	1	100	1	0.06	58	1	-	1/1	-	80	- 86	WWCNDGR
1675.872109	1675.8009	1	100	1	0.07	42	2	-	-	0/1	116	- 130	IVSDGNGMNAWVAWR
1045.614758	1045.542498	1	100	1	0.07	69	14	-	-	-	135	- 143	GTDVQAWIR

[View all couples \(1 outliers\)](#)

```

1  -----kv fgrccelaam krHGLDNYRg yslgnwvcaa kFESNFNTQA TNrNTDGSTD YGILQINSRW
81 WcNDGRtpgs rnlcnipcsa llssditasv ncaKkIVSDG NGMNAWVAWR nrckGTDVQA WIRgcr1

```

Search Engine for protein identification: MASCOT

(www.matrixscience.com)

[HOME](#)
[MASCOT](#)
[HELP](#)
[WHAT'S NEW](#)
[PRODUCTS](#)
[SUPPORT](#)
[SITE SEARCH](#)
[LINKS](#)
[EMPLOYMENT](#)
[CONTACT US](#)

Mascot: Peptide Mass Fingerprint

[Your name](#) [Email](#)

[Search title](#)

[Database](#)

[Taxonomy](#)

[Enzyme](#) [Allow up to](#) [missed cleavages](#)

[Fixed modifications](#)

[Variable modifications](#)

[Protein mass](#) [kDa](#) [Peptide tol. ±](#) [Da](#)

[Mass values](#) ☒ [MH⁺](#) ☐ [M_r](#) [Monoisotopic](#) ☒ [Average](#) ☐

[Data file](#) [Browse...](#)

[Query](#)
NB Contents of this field are ignored if a data file is specified.

[Overview](#) ☐ [Report top](#) [hits](#)

[Start Search ...](#) [Reset Form](#)

Search Engine for protein identification: ProFound

(<http://prowl.rockefeller.edu>)

ProFound
<http://prowl.rockefeller.edu/cgi-bin/ProFound>

Protein identification by comparison of a peptide map
to NCBI's nr database using a Bayesian algorithm

[Introduction](#)

[References](#)

[Example](#)

[Feedback](#)

New A new

ProFound

with more

features

is
available.

More.Digests

Set Tags

Search ProFound

Reset Form

Taxonomic category: ☐ Show species

Search for:

Approximate protein mass: - kDa

Bracket an SDS-PAGE mass by at least +/- 25%

Number of top candidate proteins:

Digest chemistry:

Note: CNBr assumes lactone product

Cysteine modified by:

Maximum number of missed cleavage sites:

Peptide masses **with mass tolerances** :

Average+/- .. Monoisotopic+/-

Average:

Monoisotopic:



Internet

PEPTIDE FINGERPRINTING SOFTWARE (1)

MW & pI: information about 2D el.

Peak list: lits of MW of peptides from in gel digestion

Sample

Sample name

Mw

pI

Peak list

Upload a file

File type

PEPTIDE FINGERPRINTING SOFTWARE (2)

Database to perform the search

Source of the protein

Include in the search only proteins in a MW range

Include in the search also splicing variants or fragments of a protein present in a database

Include in the search only proteins in a pI range

Protein

Database(s)

Taxon(s)

NCBI TaxID(s)

Keyword(s)

Swiss-Prot

TrEMBL

All

..Viruses

..Bacteria

....Escherichia coli

....Firmicutes

.....Bacillus subtilis

.....other Firmicutes

Mw min

Mw max

pI min

pI max

Varsplices

Fragments

Or

Multiple selections are possible by pressing the "Ctrl" key.

0

100000

0

14

☒

☒

PEPTIDE FINGER PRINTING SOFTWARE (3)

Number of missing cleavages allowed for each peptide

Which PTM are kept into consideration

Enzyme used for the digestion

Peptide

Digestion	
Enzyme	Trypsin
Missed cleavage	1
Resolution	Monoisotopic
Ion mode	[M+H]

PTM	
PTM	<div>Observed</div> <div>By similarity</div> <div>Probable</div> <div>Potential</div>
Max per peptide	2

Mass provided by the mass spectrometer (average or monoisotopic)

Ion modality: molecular ion [M+H] or molecular weight [M]

Modifications						
Supp	Label	Locus	Add	Supp	Mode	Max
<input type="checkbox"/>	CAM	C	C2H4ON	H	Fixed	1
<input type="checkbox"/>	MSO	M	O		Variable	2

Oxidation Trp. His (Ox(WH))
 Propionamide (PAM)
 Phosphorylation (P(ST))
 Phosphorylation (P(Y))
 Ion mode (H+)
 User defined ()

Add to list

Other modifications

Post-translational modifications allowed to the peptide

PEPTIDE FINGER PRINTING SOFTWARE (4)

Thresholds


Spectrometer			
<u>Shift max</u>	±	<input type="text" value="0.2"/>	(Dalton)
<u>Slope max</u>	±	<input type="text" value="200"/>	(ppm)
<u>Internal error max</u>	±	<input type="text" value="25"/>	(ppm)

Information about precision of the values obtained with mass spectrometer

Process Filters

<u>Minimum number of Hits</u>	<input type="text" value="4"/>
<u>pValue max</u>	<input type="text" value="1e-5"/>

Information about the identification process.
Minimum number of peptides that must be mapped on the sequence of the protein-exclusion from the results of those proteins who do not have a sufficient score



Post Translational Modifications

PTMs

Post-translational modification	Mass difference (Da)
Methylation	14.03
Propylation	42.08
Sulfation	80.06
Phosphorylation	79.98
Glycosylations by:	
Deoxyhexoses (Fuc)	146.14
Hexosamines (GlcN, GalN)	161.16
Hexoses (Glc, Gal, Man)	162.14
N-Acetylhexosamines (GlcNAc, GalNAc)	203.19
Pentoses (Xyl, Ara)	132.12
Sialic acid (NeuNAc)	291.26
Reduction of a disulfide bridge	2.02
Carbamidomethylation	57.03
Carboxymethylation	58.04
Cysteinylation	119.14
Ethylpyridylation	105.12
Acetylation	42.04
Formylation	28.01
Biotinylation	226.29
Farnesylation	204.36
Myristoylation	210.36
Pyridoxal phosphate Schiff condensation	231.14
Stearoylation	266.47
Palmitoylation	238.41
Lipoylation	188.30
Carboxylation of Asp or Glu	44.01
Deamidation of Asn or Gln	0.98
Hydroxylation	16.00
Met Oxidation	16.00
Proteolysis of a peptide bond	18.02
Deamination from Gln to pyroglutamic	-17.03



Arg, Lys



Ser, Thr, Tyr



Cys



Lys

Most common PTM

Tutte queste modificazioni post-traduzionali alterano il peso molecolare della proteina, alcune di esse ne modificano anche il punto isoelettrico: Acetilazione e Fosforilazione
 Acetilazione: rimuove una carica positiva
 Fosforilazione: introduce due cariche negative

hHMG1a

[1-106] mass = 11544.8

Cleavage at KR

Small polar:	D(1)	E(14)	N(1)	Q(6)				
Large polar:	K(16)	R(11)	H(0)					
Small non-polar:	S(14)	T(8)	A(4)	G(11)				
Large non-polar:	L(3)	I(1)	V(3)	M(0)	F(0)	Y(0)	W(0)	
Special:	C(0)	P(13)						

1 S E S S S K s s q p l a s k Q E K d g t e k R g r G R P R k 30
 31 Q P P V S P G T A L V G S Q K e p s e v p t p k R P R g r p 60
 61 k G S K n k G A A K t r K t t t t p g r K P R g r p k K l e 90
 91 k E E E E G I S Q E S S E E E Q 106

(1)	[1-6] = 623.6	(2)	[7-14] = 816.9	(3)	[15-17] = 403.4
(4)	[18-22] = 548.5	(5)	[23-23] = 174.2	(6)	[24-25] = 231.3
(7)	[26-29] = 484.6	(8)	[30-30] = 146.2	(9)	[31-45] = 1465.7
(10)	[46-54] = 983.1	(11)	[55-57] = 427.5	(12)	[58-61] = 456.5
(13)	[62-64] = 290.3	(14)	[65-66] = 260.3	(15)	[67-70] = 345.4
(16)	[71-72] = 275.3	(17)	[73-73] = 146.2	(18)	[74-80] = 732.8
(19)	[81-83] = 399.5	(20)	[84-87] = 456.5	(21)	[88-88] = 146.2
(22)	[89-91] = 388.5	(23)	[92-106] = 1738.6		

hHMG1a

[1-106] mass = 11544.8

Cleavage at KR

 UNMODIFIED
PROTEIN

Small polar:	D(1)	E(14)	N(1)	Q(6)				
Large polar:	K(16)	R(11)	H(0)					
Small non-polar:	S(14)	T(8)	A(4)	G(11)				
Large non-polar:	L(3)	I(1)	V(3)	M(0)	F(0)	Y(0)	W(0)	
Special:	C(0)	P(13)						

```

1  S E S S S K s s q p l a s k Q E K d g t e k R g r G R P R k 30
31 Q P P V S P G T A L V G S Q K e p s e v p t p k R P R g r p 60
61 k G S K n k G A A K t r K t t t t p g r K P R g r p k K l e 90
91 k E E E E G I S Q E S S E E E Q 106

```

(1)	[1-6] = 623.6	(2)	[7-14] = 816.9	(3)	[15-17] = 403.4
(4)	[18-22] = 548.5	(5)	[23-23] = 174.2	(6)	[24-25] = 231.3
(7)	[26-29] = 484.6	(8)	[30-30] = 146.2	(9)	[31-45] = 1465.7
(10)	[46-54] = 983.1	(11)	[55-57] = 427.5	(12)	[58-61] = 456.5
(13)	[62-64] = 290.3	(14)	[65-66] = 260.3	(15)	[67-70] = 345.4
(16)	[71-72] = 275.3	(17)	[73-73] = 146.2	(18)	[74-80] = 732.8
(19)	[81-83] = 399.5	(20)	[84-87] = 456.5	(21)	[88-88] = 146.2
(22)	[89-91] = 388.5	(23)	[92-106] = 1738.6		



[Open]

hHMG1a

[1-106] mass = 11784.7

Cleavage at KR

PROTEIN MODIFIED
WITH 3 PHOSPHATE
GROUPS

Small polar:	D(1)	E(14)	N(1)	Q(6)			
Large polar:	K(16)	R(11)	H(0)				
Small non-polar:	S(14)	T(8)	A(4)	G(11)			
Large non-polar:	L(3)	I(1)	V(3)	M(0)	F(0)	Y(0)	W(0)
Special:	C(0)	P(13)					

S[98] + 80.0

S[101] + 80.0

S[102] + 80.0


PHOSPHORILATED SITES: + 80 Da

1	S	E	S	S	S	K	s	s	q	p	l	a	s	k	Q	E	K	d	g	t	e	k	R	g	r	G	R	P	R	k	30
31	Q	P	P	V	S	P	G	T	A	L	V	G	S	Q	K	e	p	s	e	v	p	t	p	k	R	P	R	g	r	p	60
61	k	G	S	K	n	k	G	A	A	K	t	r	K	t	t	t	t	p	g	r	K	P	R	g	r	p	k	K	l	e	90
91	k	E	E	E	E	G	I	S	Q	E	S	S	E	E	E	Q														106	

(1)	[1-6] = 623.6	(2)	[7-14] = 816.9	(3)	[15-17] = 403.4
(4)	[18-22] = 548.5	(5)	[23-23] = 174.2	(6)	[24-25] = 231.3
(7)	[26-29] = 484.6	(8)	[30-30] = 146.2	(9)	[31-45] = 1465.7
(10)	[46-54] = 983.1	(11)	[55-57] = 427.5	(12)	[58-61] = 456.5
(13)	[62-64] = 290.3	(14)	[65-66] = 260.3	(15)	[67-70] = 345.4
(16)	[71-72] = 275.3	(17)	[73-73] = 146.2	(18)	[74-80] = 732.8
(19)	[81-83] = 399.5	(20)	[84-87] = 456.5	(21)	[88-88] = 146.2
(22)	[89-91] = 388.5	(23)	[92-106] = 1978.6		

MODIFIED PEPTIDE:
MW CHANGES

Identifying Modifications



1	AcNH ————— R	Predicted mass + 42
2	————— R	Predicted mass OK
3	————— R	Predicted mass OK
4	————— K	Predicted mass OK
5	(— K)	Probably no peptide
6	— KP ————— R	Missed cleavage
7	————— X ————— R	"Wrong" mass
8	————— K	Predicted mass OK
9	————— COOY	"Wrong" mass